

New Way of Separation of Oligonucleotides with only H₂O/MeCN in Mobile Phase

Many organic molecules of medical, pharmaceutical, industrial, and environmental interest can exist in an ionic form at specific pH levels. One such example is Oligonucleotides (OGNs), which are an important group of negatively charged molecules used in a variety of applications including medicine, research, and forensics. Liquid chromatography (LC) separation of this class of molecules is an important step in their production, characterization, and purification. There are several separation modes based on different properties of the oligonucleotide molecules.





We recently reported on a new mode of liquid chromatography which we named Bridge Ion Separation Technique (BIST[®]). Oligonucleotides are a type of molecules that can be efficiently separated by this new technique.

The method requires some amount of organic modifier in the mobile phase and a double-charged positive ion to be dynamically loaded on a cation exchange stationary phase (Fig. 1) to create the bridge formation between negatively charged stationary phase and negatively charged oligonucleotides. We found that the best result was obtained when Mg²⁺ was used as double-charged positive ions for the bridge formation.

An important property of BIST[®] mode separation is the amount of water in the MP, which affects the level of solvation of the double-charged ions and correspondingly the degree of retention of oligonucleotides. By reducing the amount of the organic modifier, any oligonucleotide can be eluted from the column Fig. 2, 3.



It turns out that the double-charged ions do not need to be supplied continuously by the MP. Once loaded with such ions the column can be operated with only $H_2O/MeCN$ mobile phase for many hundreds of injections Fig. 4. Because of the buffer free (ion-free) mobile phase, the bridging double-charged ions have very low mobility and can not be easily removed from the stationary phase, however they still provide bridge formation for the many repetitive analyses.



Copyright © 2024 SIELC Technologies, Inc. All rights reserved.





The retention of oligonucleotides in general correlates with the number of nucleobases. This is in agreement with bridge formation phenomena, since the longer chain of the molecule will have more negative charges available for bridging and longer retention will be observed. Another aspect that affects retention and the bridge formation is the geometry of the oligo-molecules. This is clearly demonstrated that the same length sequences provide very different retentions when their nucleobases are different Fig. 4, 6.



For BIST[®] to work three conditions need to be met:

- Cation exchange column with embedded negatively charge functional groups
- A double-charged positive ions should be loaded on the column
- Reduced water in the mobile phase to control the level of solvation of the bridging ions



Summary

One important advantage of this BIST[®] separation if compared to an ion-exchange mode is complete removal of the buffer from the mobile phase. When compared to a Reverse-Phase Ion-Pairing separation of oligonucleotides, the advantages are a simple MP composition, quick equilibration time with gradient mode and no interference with other molecules retained by RP mechanism. BIST[®] also offers different selectivity when compared with any other separation mode.

Every new separation mode offers new possibilities in complex analytical or prep situations.

OligoMg columns available in many standard sizes with particles 2, 3, and 5 µm and pores size 110 Å.

Copyright © 2024 by SIELC Technologies Corporation.

The information and disclosures contained herein were originated by and are the property of SIELC Technologies Corporation and may not be copied, distributed, republished, translated, uploaded, posted, transmitted, stored, or otherwise used in any shape, manner, or form by any means without the specific prior written permission of SIELC Technologies Corporation. Modification or use of this material for any reason shall be construed as a violation of the property rights of SIELC Technologies Corporation. SIELC Technologies Corporation reserves all rights, including without limitation all patent, proprietary design, manufacturing, reproduction, use, advertising, and sales rights thereto, and to any articles disclosed herein, except to the extent such rights are granted to others by law or by SIELC Technologies Corporation. BIST[®] is a registered trademark of SIELC Technologies Corporation. U.S. Patents Pending. All data were obtained in SIELC Technologies labs.